MUNDSCHUTZ TYP IIR, 3-LAGIG | PP

Produktinformation zu Art.-Nr. 291761

■ PRODUKTBESCHREIBUNG

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Produktinformation zu Art.-Nr. 291761

NELSON LABORATORIES

Bacterial Filtration Efficiency (BFE) and Differential Pressure (Delta P) Final Report

Study Number: 99918-001
Study Received Date: 27 Oct 2017
Testing Facility: Nelson Laboratories, LLC, a Business Unit of Steris Corporation, 6280 S. Redwood Rd., Salt Lake City, UT 84120 U.S.A.
Test Procedure(s): Standard Test Protocol (STP) Number: STP0102
Deviations: None

Percent BFE (%)

99.5
99.8
>99.9
99.7
99.9

Summary: The BFE test is performed to determine the filtration of the bacterial control counts upstream of the test article to the back of the Stereophococcus aureus was pre-sterilized using a nebulizer and of flow rate and test pressure. The challenge delivery was made in units (CFU) with a mean particle size (MPS) of 3.0 ± 0.5 µm. The stage, rigid particle, Accutest sample for air collection. This test EN 14683 2014, Annex B, and AS4381 2015.

The Delta P test is performed to determine the breathability of test air pressure on either side of the test article using a manometer. All was designed to comply with MIL-M-38544C, Section 4.4.1.2 and 4.4.1.2.1 and AS4381 2015.

All test method acceptance criteria were met. Testing was performed, manufacturing practice (GMP) regulations 21 CFR Parts 210, 211, and 8.

Test Side: Inside
BFE Test Area: ~40 cm²
BFE Flow Rate: 28 ± 3 L/min per minute (L/min)
Delta P Flow Rate: 8 L/min
Conditioning Parameters: 85 ± 0% relative humidity (RH) and 21 ± 0°C for a min.
Test Article Dimensions: ~170 mm x ~155 mm
Positive Control Average: 2.1 x 10⁷ CFU
Negative Monitor Count: <1 CFU
MPS: 3.0 µm

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NELSON LABORATORIES

Synthetic Blood Penetration Resistance Final Report

Study Number: 99918-001
Study Received Date: 23 Oct 2017
Testing Facility: Nelson Laboratories, LLC, a Business Unit of Steris Corporation, 6280 S. Redwood Rd., Salt Lake City, UT 84120 U.S.A.
Test Procedure(s): Standard Test Protocol (STP) Number: STP0012 Rv
Deviations: None

Synthetic Blood Penetration

None Seen

Summary: This procedure was performed to evaluate surgical face mask coating materials designed to protect against fluid penetration. The simulate an arterial spray and evaluate the effectiveness of the test possible exposure to blood and other body fluids. The distance from the cannula is 30.5 cm. A test volume of 2 mL of synthetic blood was injected.

This test method was designed to comply with ASTM F1882 EN 14683 2014 and AS4381 2015) with the following exception: performed in an environment with a temperature of 21 ± 0°C and a testing was performed at ambient conditions within one minute chamber test at those parameters.

All test method acceptance criteria were met. Testing was performed, manufacturing practice (GMP) regulations 21 CFR Parts 210, 211, and 8.

Number of Test Articles Tested: 32
Number of Test Articles Passed: 30
Test Side: Outside
Pre-Conditioning: Minimum of 4 hours at 21 ± 0°C and 85% RH
Test Conditions: 21°C and 24% RH

Results: Per ASTM F1882 and ISO 22869, an acceptable quality limit of 4.0% is met. Sampling plan when 100 of 32 test articles show passing results.

Test Pressure: 120 mmHg (16.0 kPa)
Test Article Number: Synthetic Blood Penetration
1-17-20-32
18-19
Yes

Study Director: Brandon L. Williams Study Completion Date: 05-NOV-2017

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Unterzeichnet für und im Namen von der Franz Mensch GmbH